

FUJIFILM

Wako

Code No. 165-29681 (10 μ g)
161-29683 (1 mg)

**PDGF-AA, Human, recombinant,
Animal-derived-free**

**PDGF-AA, ヒト, 組換え体,
動物由来物フリー**

(rhPlatelet-Derived Growth Factor-AA, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed human PDGF-AA

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 0.1% Trifluoroacetic Acid.

Molecular weight : 28,900 (dimer)

Biological activity : ED₅₀ is below 50 ng/mL determined by
the NR6R-3T3 cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*

This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 0.1 EU/ μ g

【Expiration date】

Indicated on the label.

【Storage】

Store at -20°C.

【Package】

Code No.	Packaging
165-29681	10 μ g
161-29683	1 mg

【Caution】

Centrifuge the vial prior to opening.

Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

Group Companies



Distributors



2505KA1