

FUJIFILM

Wako

Code No. 166-27891 (100 mg)

**Potassium Eu-encapsulated
Preyssler-type Phosphotungstate**
Eu-Preyssler 型リンタングステートカリウム塩

Potassium Eu-encapsulated Preyssler-type Phosphotungstate is a negative stain reagent. Negative staining is a qualitative method for examining the structure of proteins, bacteria, and viruses. This product with a negative staining is always obtained by a simple procedure of dissolving it in water upon use. It can be use more clearly observed compared to Potassium Na-encapsulated Preyssler-type Phosphotungstate

Appearance : White to slightly brown, crystals to powder

Solubility : Soluble in water

Molecular formula : $K_{12}[EuP_5W_{30}O_{110}] \cdot 25H_2O$

Molecular weight : 8501.53

[Storage]

Store at room temperature. Protect from the light.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA2