

FUJIFILM

Wako

< for Research Use Only >

Code No. 167-26721 (1 mg)

PA-IL Lectin, recombinant, Solution

PA-ILレクチン，組換え体，溶液

PA-IL is isolated from *Pseudomonas aeruginosa*, the gram negative bacterium and the opportunistic pathogen.

Source : *E. coli* expressed *Pseudomonas aeruginosa* lectin (PA-IL).

Appearance : PBS (pH 7.4).

Protein concentration : Indicated on the label.

Carbohydrate specificity : Gal α 1-3/4 Gal

Manufacture : Affinity purification by lactose agarose.

[Expiration date] Indicated on the label.

[Storage] Store at -20°C . After thawing, keep at $2-10^{\circ}\text{C}$.
Avoid the freeze-thaw cycles.

[Package]

Code No.	Packaging
167-26721	1 mg

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2