

FUJIFILM

Wako

Code No. 186-03691 (10 μ g)
182-03693 (1 mg)

**RANK Ligand soluble, Human,
recombinant, Animal-derived-free**
**可溶性 RANK リガンド, ヒト, 組換え体,
動物由来物フリー**
(rhsRANKL, Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed human RANK ligand soluble

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 10 mmol/L Sodium Phosphate, pH 7.5.

Molecular weight : 19,700

Biological activity : ED₅₀ is below 50 ng/mL determined by
the ability to induce bioluminescence in
engineered RANKL bioassay cells.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 1 EU/ μ g

[Expiration date]
Indicated on the label.

[Storage]
Store at -20°C.

[Package]

| Code No. | Packaging |
|-----------|------------|
| 186-03691 | 10 μ g |
| 182-03693 | 1 mg |

[Caution]
Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.fujifilm.co.jp>

Group Companies



Distributors



2505KA1