

FUJIFILM

Wako

Code No. 195-19071 (10 μ g)
191-19073 (50 μ g)

SCF Solution, Human, recombinant (Stem Cell Factor)

【Source】

Corynebacterium glutamicum expressed human Stem Cell Factor

【Formulation】

20 mmol/L HEPES, pH 7.7 (0.2 μ m filtered)

【Protein Concentration】

Indicated on the label

【Purity】

$\geq$ 95% (by SDS-PAGE)

【Biological activity】

The EC₅₀ is 1.0 ~ 20 ng/mL.

Measurement of phosphorylation assay of c-KIT using U2OS cells for reporter assay

【Endotoxin】

< 0.1 EU/ μ g

【Storage】

Store at -80°C.

Avoid repeating freeze-thaw.

【Package】

Code No.	Package
195-19071	10 μ g
191-19073	50 μ g

【Caution】

Centrifuge the vial prior to opening.

【Description】

This product is manufactured by Ajinomoto Co., Inc.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://flwkw.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2305KA2