

FUJIFILM

Wako

Code No. 195-19211 (25 μ g)
191-19213 (1 mg)

R-Spondin-1, Human, recombinant R- スポンジン -1, ヒト, 組換え体

Source : CHO cells expressed human R-Spondin-1

Formulation : Lyophilized from a 0.2 μ m filtered solution containing Phosphate Buffered Saline (PBS).

Molecular weight : approx. 26,800

Biological activity : The ED₅₀ is less than 50 ng/mL determined by the luciferase induction in HEK293 STF cells in the presence of mouse Wnt-3a.

Reconstitution : Reconstitute in sterile water to a concentration of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other aqueous buffers.

Endotoxin : Less than 0.1 EU/ μ g

[Expiration date] Indicated on the label.

[Storage] Store at -20°C .

[Package]

Code No.	Packaging
195-19211	25 μ g
191-19213	1 mg

[Caution] Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2306KA1