

FUJIFILM

Wako

Code No. 199-19351 (25 μ g)
195-19353 (1 mg)

**Sonic Hedgehog, Human,
recombinant**
ソニックヘッジホッグ, ヒト, 組換え体
(rhShh)

Source : *E. coli* expressed human sonic hedgehog

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 10 mmol/L Sodium Phosphate, pH 7.5.

Molecular weight : 20,200

Biological activity : ED₅₀ is below 5 μ g/mL determined by the
CCL-226 cell proliferation assay.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL. *Do not vortex.*
This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 0.1 EU/ μ g

[Expiration date]
Indicated on the label.

[Storage]
Store at -20°C.

[Package]

Code No.	Packaging
199-19351	25 μ g
195-19353	1 mg

[Caution]
Centrifuge the vial prior to opening.
Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

Group Companies



Distributors

