

FUJIFILM

Wako

Code No. 206-21581 (10 μ g)
202-21583 (1 mg)

**Thrombopoietin (6 × His tag),
Human, recombinant
トロンボポエチン (6 × His タグ), ヒト,
組換え体
rhTPO (6 × His tag)**

Source : CHO cells expressed human thrombopoietin (6 × His tag)

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 5 mmol/L Sodium Phosphate, 2
mmol/L L-Arginine, pH 7.5

Molecular weight : approx. 36,300

Biological activity : The ED₅₀ is less than 15 ng/mL determined
by the MO7e cell proliferation.

Reconstitution : Reconstitute in sterile water to a concentration
of 0.1 mg/mL.

Do not vortex.

This solution can then be diluted into other
aqueous buffers.

Endotoxin : Less than 1 EU/ μ g

【Expiration date】

Indicated on the label.

【Storage】

Store at -20°C.

【Package】

Code No.	Packaging
206-21581	10 μ g
202-21583	1 mg

【Caution】

Centrifuge the vial prior to opening.

Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2307KA1