

FUJIFILM

Wako

Code No. 207-21631 (20 μ g)
203-21633 (1 mg)

**Tumor Necrosis Factor- α , Mouse,
recombinant, Animal-derived-free**
腫瘍壊死因子 - α , マウス, 組換え体,
動物由来物フリー
(rmTNF- α , Animal-derived-free)

Produced using all non-animal reagents.

Source : *E. coli* expressed mouse tumor necrosis factor- α

Formulation : Lyophilized from a 0.2 μ m filtered solution
containing 10 mmol/L Sodium Phosphate, 50
mmol/L Sodium Chloride, pH 7.5.

Molecular weight : 17,400

Biological activity : ED₅₀ is below 0.1 ng/mL determined by
the cytolysis of mouse L929 cells in the
presence of Actinomycin D assay.

Reconstitution : Reconstitute in sterile water at 0.1 mg/mL. *Do
not vortex.*

This solution can then be diluted into other
aqueous buffers.

Endotoxin : Below 0.1 EU/ μ g

【Expiration date】

Indicated on the label.

【Storage】

Store at -20°C.

【Package】

Code No.	Packaging
207-21631	20 μ g
203-21633	1 mg

【Caution】

Centrifuge the vial prior to opening.

Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2402KA1