

FUJIFILM

Wako

Code No. 220-02041 (500 μ g)

**Vitronectin (20-398 aa), Human,
recombinant, Solution**

ビトロネクチン, ヒト, 組換え体, 溶液
(rhVTN (20-398 aa) Solution)

Source : *E. coli* expressed human vitronectin (20-398 aa)

Formulation : Sterile-filtered, in 20 mM pH 8.0 Tris-HCl
Buffer, with a proprietary formulation of NaCl,
KCl, EDTA, Arginine, DTT and Glycerol.

Protein : Indicated on the label.

Assay (SDS-PAGE) : min. 90%

[Expiration date] Indicated on the label.

[Storage] Store at -80°C .

[Package] 500 μ g.

[Caution] Avoid repeated freeze-thaw cycle.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2308KA2