

FUJIFILM

Wako

Code No. 222-02481 (10 μ g)
228-02483 (50 μ g)

VEGF-165 Solution, Human, recombinant (Vascular Endothelial Growth Factor-165)

[Source]

Corynebacterium glutamicum expressed human vascular endothelial growth factor-165

[Formulation]

20 mmol/L GSSG and 20 mmol/L sodium citrate, pH 3.5 (0.2 μ m filtered)

[Protein Concentration]

Indicated on the label

[Purity]

$\geq$ 95% (by SDS-PAGE)

[Biological Activity]

The EC₅₀ is 1 ~ 30 ng/mL

Measurement of VEGFR2/NFAT reporter assay

[Endotoxin]

< 0.1 EU/ μ g

[Storage]

Store at -80°C

Avoid repeating freeze-thaw.

[Package]

Code No.	Package
222-02481	10 μ g
228-02483	50 μ g

[Caution]

Centrifuge the vial prior to opening.

[Description]

This product is manufactured by Ajinomoto Co., Inc.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2310KA1