




























for Pharmaceutical, Food, and Water Environment Analysis Analytical Standards of Nitrosamines

Nitrosamines are compounds of the chemical structure where hydrogen of the amine nitrogen is substituted with a nitroso group. And these are known to be carcinogenic. Nitrosamine are often detected as impurities during the manufacturing process of a drug product. Recently, nitrosamines were detected in some sartan and ranitidine drug products, and these products were recalled. In response to this event, in September 2019, the EMA (European Medicines Agency) requested the marketing authorization holders to evaluate the risk of the presence of nitrosamine impurities in their drugs and take appropriate risk mitigation measures.



Product List

<i>N</i> -Nitrosodimethylamine Standard 		<i>N</i> -Nitrosomethylethylamine Standard 		<i>N</i> -Nitrosodiethylamine Standard 	
for Gas Chromatography		for Chromatography		for Chromatography	
					
Code No. 147-03781	Vol. 1 g	Code No. 140-10001	Vol. 50 mg	Code No. 141-09921	Vol. 100 mg
Abbrev. NDMA	CAS RN® 62-75-9	Abbrev. NEMA, NMEA	CAS RN® 10595-95-6	Abbrev. NDEA	CAS RN® 55-18-5
<i>N</i> -Nitrosoethylisopropylamine Standard 		<i>N</i> -Nitrosodi- <i>n</i> -propylamine Standard 		<i>N</i> -Nitrosodiisopropylamine Standard 	
for Chromatography		for Chromatography		for Chromatography	
					
Code No. 142-09951	Vol. 50 mg	Code No. 140-09991	Vol. 100 mg	Code No. 145-09941	Vol. 50 mg
Abbrev. EIPNA, NEIPA	CAS RN® 16339-04-1	Abbrev. NDPA	CAS RN® 621-64-7	Abbrev. DIPNA, NDIPA	CAS RN® 601-77-4
<i>N</i> -Nitrosodi- <i>n</i> -butylamine Standard 		<i>N</i> -Nitrosomethylaminobutyric Acid Standard 		<i>N</i> -Nitrosomethylphenylamine Standard 	
for Chromatography		for Chromatography		for Chromatography	
					
Code No. 149-09961	Vol. 100 mg	Code No. 146-09971	Vol. 50 mg	Code No. 148-09931	Vol. 50 mg
Abbrev. NDBA	CAS RN® 924-16-3	Abbrev. NMBA	CAS RN® 61445-55-4	Abbrev. NMPA	CAS RN® 614-00-6
<i>N</i> -Nitroso- <i>N'</i> -methylpiperazine Standard 		<i>N</i> -Nitrosomorpholine Standard 		<i>N</i> -Nitrosodimethylamine- <i>d</i> ₆ Standard 	
for Chromatography		for Chromatography		for Chromatography	
					
Code No. 143-09981	Vol. 50 mg	Code No. 141-10031	Vol. 100 mg	Code No. 144-10021	Vol. 100 mg
Abbrev. MeNP	CAS RN® 16339-07-4	Abbrev. NMOR	CAS RN® 59-89-2	Abbrev. NDMA- <i>d</i> ₆	CAS RN® 17829-05-9
<i>N</i> -Nitrosodiethylamine- <i>d</i> ₁₀ Standard 					
for Chromatography					
					
Code No. 147-10011	Vol. 50 mg				
Abbrev. NDEA- <i>d</i> ₁₀	CAS RN® 1219794-54-3				

Regulatory Trends in Each Country

Nitrosoamine	Abbreviation	FDA	USP	EMA	Ph. Eur.
<i>N</i> -Nitrosodi- <i>n</i> -butylamine	NDBA	✓	✓	✓	✓
<i>N</i> -Nitrosodiethylamine	NDEA	✓	✓	✓	✓
<i>N</i> -Nitrosodiisopropylamine	DIPNA, NDIPA	✓	✓	✓	✓
<i>N</i> -Nitrosodimethylamine	NDMA	✓	✓	✓	✓
<i>N</i> -Nitrosodi- <i>n</i> -propylamine	NDPA				✓
<i>N</i> -Nitrosoethylisopropylamine	EIPNA, NEIPA	✓	✓	✓	✓
<i>N</i> -Nitrosomethylaminobutyric Acid	NMBA	✓	✓	✓	✓
<i>N</i> -Nitrosomethylphenylamine	NMPA		✓	✓	
<i>N</i> -Nitroso- <i>N'</i> -methylpiperazine	MeNP			✓	
<i>N</i> -Nitrosomorpholine	NMOR			✓	

Analytical Methods of USP and Ph. Eur

■ Ph. Eur. 10.6 2.4.42 *N*-Nitrosamine impurities in active substance (As of June 2021)

	Analyte	Internal Standard
Procedure A	LC/MS/MS	NDMA, NDEA, NMBA, DIPNA, EIPNA
Procedure B	GC/MS	NDMA, NDEA
Procedure C	GC/MS/MS	NDMA, NDEA, NDBA, DIPNA, EIPNA, NDPA

■ USP <1469> *N*-Nitrosamine Impurities (As of June 2021)

	Analyte	Internal Standard
Procedure 1	LC/HRMS	NDMA, NDEA, NDIPA, NEIPA, NMBA, NDBA, NMPA
Procedure 2	GC/MS	NDMA, NDEA, NDIPA, NEIPA
Procedure 3	LC/MS/MS	NDMA, NDEA, NDIPA, NEIPA, NMBA, NDBA
Procedure 4	GC/MS/MS	NDMA, NDEA, NDIPA, NEIPA, NDBA, NMPA

*The regulated components and analytical methods are based on the information available as of June 2021.

For details, see the website of each pharmacopoeia.

*Please check our website for details.

Pharma Manufacturing & QC > QC Test > Analytical Standards of Nitrosamines

<https://labchem-wako.fujifilm.com/us/category/02128.html>

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: <https://labchem-wako.fujifilm.com/> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Laboratory Chemicals site

<https://labchem-wako.fujifilm.com>



FUJIFILM Wako Pure Chemical Corporation
1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan
Tel: +81 6 6203 3741 Fax: +81 6 6203 1999
ffwk-cs@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation
1600 Bellwood Road, Richmond, VA 23237, U.S.A.
Toll-Free (U.S. only): +1 877 714 1920
Tel: +1 804 271 7677 Fax: +1 804 271 7791
wkuslabchem@fujifilm.com

FUJIFILM Wako Chemicals (Hong Kong) Limited
Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road,
Tsuen Wan, N.T., Hong Kong
Tel: +852-2799-9019 Fax: +852-2799-9808
wkhk.info@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH
Fuggerstr 12, 41468 Neuss, Germany
Tel: +49 2131 311 0 Fax: +49 2131 311 100
labchem_wkeu@fujifilm.com

FUJIFILM Wako (Guangzhou) Trading Corporation
Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road,
Guangzhou, 510095, China
Tel: +86-20-8732-6381(Guangzhou) Tel: +86-21-6288-4751(Shanghai)
Tel: +86-10-6413-6388(Beijing)
wkgz.info@fujifilm.com