

FUJIFILM

Wako

Code No. 016-20001 (50 µg)

Anti Iba1, Rabbit (for Western Blotting)
抗 Iba1, ウサギ (ウエスタンブロッティング用)

Rabbit Anti Iba1 antibody is raised against synthetic peptide corresponding to C-terminus of Iba1.

For research use only. Not for use in diagnostic procedure. ; use in any such application is the responsibility of the user.

Preparation : Purified by the antigen affinity chromatography from rabbit antisera and prepared in TBS solution. Contains no preservatives and stabilizers.

Specificity : Reactive with human, mouse and rat Iba1.

Working Concentration : Westernblot 0.5-1 µg/mL

Storage : Keep at - 20°C.

After opening, aliquot contents and freeze at - 20°C.

Package : 50 µg (100 µL)

References :

- 1) Imai, Y., Ibata, I., Ito, D., Ohsawa, K. and Kohsaka, S. : *Biochem. Biophys. Res. Commun.*, **224**, 855 (1996).
- 2) Ito, D., Imai, Y., Ohsawa, K., Nakajima, K., Fukuuchi, Y. and Kohsaka, S. : *Brain Res. Mol. Brain Res.*, **57**, 1 (1998).
- 3) Ohsawa, K., Imai, Y., Kanazawa, H., Sasaki, Y. and Kohsaka, S. : *J. Cell Sci.*, **113**, 3073 (2000).
- 4) Sasaki, Y., Ohsawa, K., Kanazawa, H., Kohsaka, S. and Imai, Y. : *Biochem. Biophys. Res. Commun.*, **286**, 292 (2001).
- 5) Kanazawa, H., Ohsawa, K., Sasaki, Y., Kohsaka, S. and Imai, Y. : *J. Biol. Chem.*, **277**, 20026 (2002).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-311-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2003KA4