
New Product

GREEN CHEMISTRY

Supplement toWako Product Update No. 9 —

Contents

A. Alphabetical Index	2
B. Metal Catalyst Index	4
C. New Products	5
1. Metal Catalysts	5
a. OsO ₄ , Immobilized Catalyst I	5
b. Polymer-Incarcerated Pd	6
c. 3D-Vacant Pd Complexes	7
d. Chiral Zr Catalyst with Molecular Sieves	8
e. Precious Powder Metal Catalysts	9
f. Metal Scavengers	9
2. lonic Liquid	10
3. Chiral Brønsted Acid Catalyst	11
4. Rare Earth Triflate	11
5. Maruoka Catalysts	12



Please visit the Wako Online Catalog http://search.wako-chem.com



Green Chemistry A. ALPHABETICAL INDEX

	WPU#-page	Description	Category	Metal
Α	9-27	Acetone, Dehydrated	Dehydrated Solvent	-
	9-27	Acetonitrile, Dehydrated	Dehydrated Solvent	-
	9-1,9	Acetonitrile Dichloropalladium (II)-Triphenylphosphine PE fibres	Polymer-supported Catalyst	Pd
	9-20	Acylase, Amano	Biocatalyst	-
	9-12	linear-Alkylbenzenesulfonic Acid	Reaction in Water	-
	9-20	D-Aminoacylase Amano	Biocatalyst	-
	9-21	Anthracene-9,10-bis(5-resorcinol)	Anthracene	-
В	9-27	Benzene, Dehydrated	Dehydrated Solvent	-
	9-22	(R)-(+)-BINAP, (S)-(-)-BINAP	Asymmetric Ligand	-
	9-22	(R)-2,2'-Binaphthyl-14-crown-4, (S)-2,2'-Binaphthyl-14-crown-4, etc.	Asymmetric Ligand	-
	9-23	(R,R)-linked-BINOL and (S,S)-linked-BINOL	Asymmetric Ligand	_
	9-1,9	cis-Bis(acetonitrile)dichloroplatinum (II)-Triphenylphosphine PE fibres	Polymer-supported Catalyst	Pt
	9-22			-
		(R)-(+)-2,2'-Bis (diphenylphosphino)-1,1'-binaphthyl, etc.	Asymmetric Ligand	
	14-12	(R,R)-3,5-Bisfluoromethylphenyl-NAS Bromide	Chiral Phase-Transfer Catalysts/Maruoka Catalysts	-
	14-11	(R)-(-)-3,3'-Bis(4-nitrophenyl)-1,1'-binaphthyl-2,2' diyl Hydrogen Phosphate	Non-metal Catalyst/Chiral BrΦnsted Acid Catalyst	-
-	9-2,9	Bis(η-norbornadiene)rhodium (I) Tetrafluoroborate-Triphenylphosphine PE fibres	Polymer-supported Catalyst	Rh
	9-13	(R,R)-3,5-Bistrifluoromethylphenyl-NAS Bromide	Chiral Phase-Transfer Catalyst	-
	9-27	1-Butanol, Dehydrated	Dehydrated Solvent	-
	9-27	2-Butanone, Dehydrated	Dehydrated Solvent	-
	9-27	Butyl Acetate, Dehydrated	Dehydrated Solvent	-
C	9-27	Chloroform, Dehydrated	Dehydrated Solvent	-
	9-30	Chromato Sheet	Analytical Chemistry	
	9-1,16	Cobaltocene	Metallocene/Catalyst	Co
	9-27	Cyclohexane, Dehydrated	Dehydrated Solvent	-
D	9-12	DBSA	Reaction in Water	-
_	9-1,4	Dibromobis(triphenylphosphine)nickel (II), Supported PS Resin	Polymer-supported Catalyst	Ni
	9-1,5	Di-µ-chlorobis [(η-allyl) palladium (II)], Supported PEG-PS Resin	Polymer-supported Catalyst	Pd
	9-2,9	Di-μ-chlorobis (η-norbonadiene)dirhodium (I)-Triphenylphosphine PE fibres	Polymer-supported Catalyst	Rh
	9-1,4	Dichlorobis(triphenylphosphine)cobalt(II), Supported PS Resin	Polymer-supported Catalyst	Со
-	9-27	Dichloromethane, Dehydrated	Dehydrated Solvent	-
	9-27	Diethyl Ether, Dehydrated	Dehydrated Solvent	-
	9-1,9	Dihydrogen Hexachloroplatinate (IV) n-Hydrate-Pyridine PE fibres	Polymer-supported Catalyst	Pt
	9-25	4-(4,6-Dimethoxy-1,3,5-triazin-2-yl)-4-methylmorpholinium Chloride n-Hydrate	Condensing agent in Water	-
	9-27	N,N-Dimethylacetamide, Dehydrated	Dehydrated Solvent	-
	9-27	N,N-Dimethylformamide, Dehydrated	Dehydrated Solvent	-
	9-27	Dimethyl Sulfoxide, Dehydrated	Dehydrated Solvent	-
	9-27	1,4-Dioxane, Dehydrated	Dehydrated Solvent	-
	9-1,9	Dipotassium Dioxotetrahydroxoosmate (VIII)-Triethylamine PE fibres	Polymer-supported Catalyst	Os
	9-2,8	met-DIRUX	Chiral Catalyst	Ru
	9-25	DMT-MM	Condensing agent in Water	
	9-12	Dodecylbenzenesulfonic acid	Reaction in Water	-
	9-24	Dodecyldimethylsulfonium lodide	Alternative to DMSO	+
		Dodecyl Methyl Sulfide		_
	9-24		Alternative to dimethylsulfide	
_	9-24	Dodecyl Methyl Sulfoxide	Alternative to dimethylsulfoxide	-
E	9-27	Ethanol, Dehydrated	Dehydrated Solvent	-
	9-27	Ethyl Acetate, Dehydrated	Dehydrated Solvent	-
	14-7	(Ethylenediamine)dinitratopalladium(II) <[en)Pd](ONO ₂) ₂]>	3D water-soluble complex	Pd
	9-27	Ethylene Glycol, Dehydrated	Dehydrated Solvent	-
	9-21	1-Ethyl-3-methylimidazolium Trifluoromethanesulfonate	lonic Liquid	-
F	9,1,9	FibreCat®-Os	Immobilized Precious Metal Catalyst	Os
	9-1,9	FibreCat®-Pd	Polymer-supported Catalyst	Pd
	9-1,9	FibreCat®-Pt	Immobilized Precious Metal Catalyst	Pt
	9-1,9	FibreCat®-Rh	Immobilized Precious Metal Catalyst	Rh
	9-1,9	FibreCat®-Ru	Immobilized Precious Metal Catalyst	Ru
н	9-p.26	HCFC-225 (mixture of CF ₃ CF ₂ CHCl ₂ and CCIF ₂ CF ₂ CHClF)	CFC-Alternatives/Solvent	-
	9-27	Heptane, Dehydrated	Dehydrated Solvent	-
	9-27	Hexane, Dehydrated	Dehydrated Solvent	-
ı	14-8	(R)-3-I-ZrMS [MS: molecular sieves]	Molecular Sieves-supported Catalyst	Zr
	14-11	Lanthanam(III) Trifluoromethanesulfonate	Reaction in water/Rare-earth Triflate	La
L	9-23	(R,R)-linked-BINOL, (S,S)-linked-BINOL	Asymmetric Ligand	_

Green Chemistry A. ALPHABETICAL INDEX

	WPU#-page	Description	Category	Metal
М	9-13; 14-12	Maruoka Catalyst	Chiral Phase-Transfer Catalyst	-
	9-1 ,16	Metallocene-Cr	Catalyst/Metallocene	Cr
	9-1 ,16	Metallocene-Hf	Catalyst/Metallocene	Hf
	9-2 ,16	Metallocene-Ti	Catalyst/Metallocene	Ti
	9-2 ,16	Metallocene-W	Catalyst/Metallocene	w
	9-2 ,16, 17	Metallocene-Zr	Catalyst/Metallocene	Zr
	14-9	Metal Scavenger, Degussa type Deloxan®	Metal Scavenger (organofunctional polysiloxanes)	-
	9-2,8	met-DIRUX	Chiral Catalyst	Ru
	9-2,8	Methanethiolate-Bridged Diruthenium Complex	Chiral Catalyst	Ru
	9-27	Methanol, Dehydrated	Dehydrated Solvent	-
	14-10	1-(2-Methyl-18-crown-6)-3-methylimidazolium Hexafluorophosphate	Ionic Liquid	-
	9-27	4-Methyl-2-pentanone, Dehydrated	Dehydrated Solvent	-
	9-27	1-Methyl-2-pyrrolidone, Dehydrated	Dehydrated Solvent	-
N	9-18	Nafion® Dispersion Solution	Solid Super-Acid Catalyst	-
	9-18	Nafion® NR-50	Solid Super-Acid Catalyst	-
	14-11	Neodymium Trifluoromethanesulfonate	Reaction in water/Rare-earth Triflate	Nd
	14-11	(R)-NP-NAP	Non-metal Catalyst/Chiral BrΦnsted Acid Catalyst	
0	14-5	Os IC-I	Polymer-supported Catalyst	Os
	14-5	Osmium Oxide, Immobilized Catalyst I	Polymer-supported Catalyst	Os
	9-1,3	Osmium (VIII) Oxide, Microencapsulated	Polymer-supported Catalyst	Os
	9-1,9	Osmium (VIII) Oxide-Pyridine PE fibres	Polymer-supported Catalyst	Os
	9-23	3,3"-[Oxybis(methylene)]bis-(1R,1"R)-1,1'-bi-2-naphthol, etc.	Asymmetric Ligand	-
P	9-1,9	Palladium (II) Acetate-Dicyclohexylphenylphosphine PE fibres	Polymer-supported Catalyst	Pd
-	9-1,9	Palladium (II) Acetate-Dicyclonexyphenyiphosphine PE fibres	Polymer-supported Catalyst	Pd
		Palladium-Activated Carbon Ethylenediamine Complex	Stable Pd/C Catalyst	Pd
	9-1,6		·	Pd
	9-1,5	Palladium (II)-Hydrotalcite	Polymer-supported Catalyst	
	14-7	Palladium-Nanobowl	3D water-soluble complex	Pd
	9-1 ,6, 11; 14-7	Palladium-Nanocage	Catalyst acted in Water	Pd
	14-9	Palladium on Activated Carbon, Degussa	Precious Metal Powder Catalysts	Pd
	9-1,4	PEP-Pd	Polymer-supported Catalyst	Pd
	9-1,4	PEG-PS resin-supported phosphine-Palladium complex	Polymer-supported Catalyst	Pd
	14-6	PI Pd (Pd abt. 3%)	Polymer-incarcerated Catalyst	Pd
	9-26	Poly[p-(hydroxy)(tosyloxy)iodostyrene]	Polymer-supported reagent	-
	14-6	Polymer-Incarcerated Palladium	Polymer-incarcerated Catalyst	Pd
	9-27	Propanol, Dehydrated	Dehydrated Solvent	-
	9-27	Pyridine, Dehydrated	Dehydrated Solvent	-
R	14-9	Rhodium on Activated Carbon, Degussa	Precious Metal Powder Catalysts	Rh
	9-2,7	RuHAP	Polymer-supported Catalyst	Ru
	9-2,7	Ruthenium(III)-Hydroxyapatite	Polymer-supported Catalyst	Ru
	9-2 ,16	Ruthenocene	Metallocene/Catalyst	Ru
S	9-2 ,12	Scandium tridodecylsulfonate	Catalyst/Reaction in water/Rare-earth Triflate	Sc
	91,4	Scandium Trifluoromethanesulfonate, Microencapsulated	Polymer-supported Catalyst	Sc
	9-1 ,11; 14-11	Scandium (III) Trifluoromethanesulfonate	Catalyst/Reaction in water/Rare-earth Triflate	Sc
	9-1,9	Sodium Ruthenate (VII)-Triethylamine PE fibres	Polymer-supported Catalyst	Ru
	9-2 ,12	STDS	Reaction in Water	Sc
Т	9-14	TaDiAS-[(4R,5R)-2-t-butyl-2-methyl-N,N,N'N'-tetrakis(4-methoxylbenzyl)] Diiodide, etc.	Chiral Phase-Transfer Catalyst	-
	9-14	TaDiAS-[(4R,5R)-2,2-dipropyl-N,N,N',N'-tetrakis(4-methylbenzyl)] Bis(tetrafluoroborate), etc.	Chiral Phase-Transfer Catalyst	-
	9-27	Tetrahydrofuran, Dehydrated	Dehydrated Solvent	-
	9-1, 4	Tetrakis(triphenylphosphine)palladium(0), Supported PS Resin	Polymer-supported Catalyst	Pd
	14-11	Thulium Trifluoromethanesulfonate	Reaction in water/Rare-earth Triflate	Tm
	9-27	Toluene, Dehydrated	Dehydrated Solvent	-
	9-13; 14-7	(R,R)-3,4,5-Trifluorophenyl-NAS Bromide	Chiral Phase-Transfer Catalysts/Maruoka Catalysts	-
	14-12	(S,S)-3,4,5-Trifluorophenyl-NAS Bromide	Chiral Phase-Transfer Catalysts/Maruoka Catalysts	-
	9-24	p-(Trimethylsilyl)benzenthiol	Less Odor Sulfur Compound	-
	9-24	p-(Trimethylsilyl)phenylmethanethiol	Less Odor Sulfur Compound	-
K	9-27	Xylene, Dehydrated	Dehydrated Solvent	-
Υ	14-11	Ytterbium(III) Trifluoromethanesulfonate n-Hydrate	Reaction in water/Rare-earth Triflate	Yb
z	9-2 ,19	Zirconia, Sulfated	Solid Super-Acid Catalyst	Zr

Please visit "http://www.wako-chem.co.jp/siyaku/journal/proup/pdf/proup9.pdf/" to browse the Wako Product Update No. 9

Green Chemistry B. Metal Catalysts Index

	Category	Description	WPU # - pag
Co	Metallocene		9-16
	Polymer-supported Catalyst	Dichlorobis(triphenylphosphine)cobalt(II), Supported PS Resin	9-4
Cr	Metallocene		9-16
Hf	Metallocene		9-16
La	Reaction in water/Rare-earth Triflate	Lanthanam(III) Trifluoromethanesulfonate	14-11
Nd	Reaction in water/Rare-earth Triflate	Neodymium Trifluoromethanesulfonate	14-11
Ni	Polymer-supported Catalyst	Dibromobis(triphenylphosphine)nickel (II), Supported PS Resin	9-4
Os	Polymer-supported Catalyst	Dipotassium Dioxotetrahydroxoosmate (VIII)-Triethylamine PE fibres	9-9
	Polymer-supported Catalyst	FibreCat® 3003, 3004	9-9
	Polymer-supported Catalyst	Osmium Oxide, Immobilized Catalyst I <osic-i></osic-i>	14-5
	Polymer-supported Catalyst	Osmium (VIII) Oxide, Microencapsulated; PEM Microencapsulated	9-3
	Polymer-supported Catalyst	Osmium (VIII) Oxide-Pyridine PE fibres	9-9
Pd	Polymer-supported Catalyst	Acetonitrile Dichloropalladium (II)-Triphenylphosphine PE fibres	9-9
	Polymer-supported Catalyst	Di-μ-chlorobis [(η-allyl) palladium (II)], Supported PEG-PS Resin	9-6
	3D water-soluble complex	(Ethylenediamine)dinitratopalladium(II) <[en)Pd](ONO ₂) ₂]>	14-7
	·	FibreCat* 1001, 1007, 1026	9-9
	Polymer-supported Catalyst		
	Polymer-supported Catalyst	Palladium (II) Acetate-Dicyclohexylphenylphosphine PE fibres	9-9
	Polymer-supported Catalyst	Palladium (II) Acetate-Triphenylphosphine PE fibres	9-9
	Polymer-supported Catalyst	Palladium (II)-Hydrotalcite	9-5
	Stable Pd/C Catalyst	Palladium-Activated Carbon Ethylenediamine Complex	9-6
	3D water-soluble complex	Palladium-Nanobowl	14-7
	3D water-soluble complex	Palladium-Nanocage	9-11, 14-
	Precious Metal Powder Catalysts	Palladium on Activated Carbon, Degussa	14-9
	Polymer-supported Catalyst	PEG-PS resin-supported phosphine-Palladium complex	9-5
	Polymer-supported Catalyst	PEP-Pd	9-5
	Polymer-incarcerated Catalyst	PI Pd	14-6
	Polymer-supported Catalyst	Tetrakis(triphenylphosphine)palladium(0), Supported PS Resin	9-4
	Polymer-supported Catalyst	Triphenylphosphine Pd (0), Microencapsulated	9-4
Pt	Polymer-supported Catalyst	cis-Bis(acetonitrile)dichloroplatinum (II)-Triphenylphosphine PE fibres	9-9
	Polymer-supported Catalyst	Dihydrogen Hexachloroplatinate (IV) n-Hydrate-Pyridine PE fibres	9-9
	Polymer-supported Catalyst	FibreCat® 4001, 4003	9-9
	Precious Metal Powder Catalysts	Platinum on Activated Carbon, Degussa	14-9
Rh	Precious Metal Powder Catalysts	Rhodium on Activated Carbon, Degussa	14-9
	Polymer-supported Catalyst	$Bis(\eta-norbornadiene) rhodium \ (I) \ Tetrafluoroborate-Triphenylphosphine \ PE \ fibres$	9-9
	Polymer-supported Catalyst	Di-μ-chlorobis (η-norbonadiene)dirhodium (I)-Triphenylphosphine PE fibres	9-9
	Polymer-supported Catalyst	FibreCat® 2003, 2006	9-9
Ru	Polymer-supported Catalyst	FibreCat® 3002	9-9
	Catalyst	met-DIRUX	9-8
	Catalyst	Methanethiolate-Bridged Diruthenium Complex	9-8
	Polymer-supported Catalyst	RuHAP	9-7
	Polymer-supported Catalyst	Ruthenium(III)-Hydroxyapatite	9-7
	Metallocene		9-16
	Polymer-supported Catalyst	Sodium Ruthenate (VII)-Triethylamine PE fibres	9-9
Sc	Reaction in Water/Rare-earth Triflate	Scandium (III) Trifluoromethanesulfonate	9-11, 14-
-	Reaction in Water	Scandium tridodecylsulfonate	9-12
	Polymer-supported Catalyst	Scandium Trifluoromethanesulfonate, Microencapsulated	9-4
	Reaction in Water	STDS	9-12
Ti	Metallocene		9-16
Tm	Reaction in water/Rare-earth Triflate	Thulium Trifluoromethanesulfonate	14-11
w	Metallocene	The state of the s	9-16
VV		Zircania Tungstata	
	Solid Super-Acid Catalyst	Zirconia, Tungstate	9-19 14-11
VL	Deposition in const. (D	Ytterbium(III) Trifluoromethanesulfonate n-Hydrate	1/1-11
Yb	Reaction in water/Rare-earth Triflate	receivant, in macrometrate and in the control of th	
Yb Zr	Metallocene		9-16, 17
		(R)-3-I-ZrMS Zirconia, Sulfated	9-16, 17 14-8 9-19

a. OsO₄, Immobilized Catalyst I [Os IC-I]

Osmium Oxide Adsorbed on Polymer Surfaces with High Solvent Resistance

Osmium (VIII) Oxide is one of the best reagents for converting olefin into diol and is in demand for various purposes. However, there are almost no successful examples of it in terms of actual industrial processes. This is due to the fact that Osmium (VIII) Oxide is highly toxic, expensive and volatile; handling of the catalyst is difficult. This results in the need for special precautions in the workplace. There are also concerns for its impact and burden on the global environment. Os IC-I features an inexpensive fixed catalyst where the volatility and toxicity are reduced by having the Osmium (VIII) Oxide absorbed on polymer surfaces. In addition, using a polymer with high solvent resistance makes it usable with various reaction solvents, allowing for easier setting of reaction conditions.

[Features]

- 1. Readily recoverable and reusable by filtration
- 2. Easy separation from reaction mixtures
- 3. Reduced toxicity and irritating odor by suppressing volatility
- 4. Various reaction solvents can be used due to high solvent resistance

[Reaction 1]

Diol formation from styrene (reoxidant: NMO)

<Reaction condition>

0.5 g (4.8 mmol) of styrene, 5 mL of CH₃CN, 5 mL of acetone, and 5 mL of water are combined, to which 0.97 g (7. 2mmol) of NMO and **Os IC-I** (OsO₄: equivalent of 0.24 mmol) are added. The mixture is stirred for 16 hours at room temperature. After completion of the reaction, while washing with ethyl acetate and water, **Os IC-I** is filtered off, the mother liquor is extracted twice with ethyl acetate, and with the organic layer the mixture is washed with 10 % sodium thiosulfate solution followed by saturated saline. In addition, after the organic layer is dried with anhydrous magnesium sulfate, filtered, and vacuum concentrated, the crude product is refined by silica gel column chromatography.

[Reaction 2]

Diol formation from styrene (reoxidant: K₃Fe(CN)₆)

Run	1	2	3
Yield of Product (%)	85	85	83
ee (%)	XX	XX	XX

<Reaction condition>

0.5 g (4.8 mmol) of styrene, 8 mL of t-BuOH, and 8 mL of water are combined, to which 1.99 g (14.4 mmol) of K_2CO_3 , 4.74 g (14.4 mmol) of $K_3Fe(CN)_6$, 187 mg (0.24 mmol) of $(DHQD)_2PHAL$, and **Os IC-I** (OsO₄: equivalent of 0.24 mmol) are added. The mixture is stirred for 16 hours at room temperature. After completion of the reaction, while washing with ethyl acetate and water, **Os IC-I** is filtered off, the mother liquor is extracted twice with ethyl acetate, and with the organic layer the mixture is washed with 10 % sodium thiosulfate solution followed by saturated saline. In addition, after the organic layer is dried with anhydrous magnesium sulfate, filtered, and vacuum concentrated, the crude product is refined by silica gel column chromatography.

[Comparison of recovery rates according to solvents used]

entry	Solv.	Yield (%)
1	water/THF(1/1)	85
2	water/CHCl ₃ (1/1)	73
3	water/AcOEt(1/1)	81
4	water/toluene/CH ₃ CN(1/1/1)	88

Note: Depending on the reaction substrates and reaction conditions, osmium may be dissolved out.


Description	abbreviation	Grade	Wako Cat No.	Package Size
Osmium Oxide, Immobilized Catalyst I	Os IC-I	for Organic Synthesis	153-02581	5 g
Osmium Oxide, immobilized Catalyst i		Tor Organic synthesis	151-02582	25 g



- Other services including industrial supply of osmium oxide, immobilized catalyst I, as custom synthesis of well as diol formation are available.
- Please free to contact us to inquire about other immobilized osmium oxides.

b. Polymer-Incarcerated Palladium [PI Pd]

 ${f PI}$ ${f Pd}$ is a new type of immobilized catalyst with the organic polymer physically supporting the palladium metal. The palladium catalyst features excellent resistance to solvents, high pressure and high temperature. Excellent activity of ${f PI}$ ${f Pd}$ has been demonstrated in hydrogenation of various olefins, benzyl ethers, and nitro and aromatic compounds. In addition, it can be recovered by simple filtration after the reaction and reused.



[Reactions]

1) Suzuki-Miyaura reaction²⁾

2) Allylic substitution reaction¹⁾

3) Hydrogenation of Benzalacetone^{1), 3)}

Description	Purpose	Wako Cat No.	Package Size
Pl Pd (Pd: abt. 3 %)	for coupling reactions	168-21991	1 g
PIPU (Pu: abt. 5 %)	for coupling reactions	164-21993	5 g

[References]


- 1) R. Akiyama, S. Kobayashi: J. Am. Chem. Soc., 125, 3412(2003).
- 2) K. Okamoto, R. Akiyama, S. Kobayashi: Org. Lett., 6, 1987(2004).
- 3) K. Okamoto, R. Akiyama, S. Kobayashi: J. Org. Chem., 69, 2871(2004).

c. 3D-Vacant Pd Complexes

Water-soluble palladium complex with 3D pore

Certain three-dimensional chemical compounds have closed spaces, created by their skeletons, which are isolated from the outside. When a molecule is incorporated into such a space, it is expected to express new physical properties and reactivity. Pd-Nanocage and Pd-Nanobowl are types of 3D complexes with nanosized pores, which are instantaneously and quantitatively constructed by self-assembly driven by coordinate bonds with a transition metal. These 3D pores can provide a space for various reactions, including specific reactions using the shape and size of the isolated space, specific clathration of organic compounds, and organic synthesis in water. In addition to Pd-Nanocage, two new products have been added to this category; Pd-Nanobowl and [(en)Pd](ONO₂)₂, which aids in the assembly of 3D complexes.

[Structural formula]




[Reactions]

Diels-Alder reaction

Wacker reaction

[2+2] Photocyclization reaction between different molecules



Description	Purpose	Wako Cat No	Package Size
Palladium-Nanocage	For organic synthesis	160-20471	1 g
Palladium-Nanobowl	For organic synthesis	165-21761	200 mg
	For organic synthesis	161-21763	1 g
(Ethylenediamine)dinitratopalladium(II)	Fau augustia armatha air	052-07341	200 mg
Abbr.: [(en)Pd] (ONO ₂) ₂	For organic synthesis	058-07343	1 g



- 1) D. Oguro, M. Miyazawa, H. Oka, k. Yamaguchi, K. Ogura, M. Fujita: Nature, 378, 469 (1995).
- 2) T. Kusukawa, M. Fujita: Angew. Chem., Int. Ed. Engl., 37, 3142 (1998).
- 3) H. Ito, T. Kusukawa, M. Fujita: Chem. Lett., 598 (2000).
- 4) M. Yoshizawa, Y. Takeyama, T. Okano, M. Fujita: J. Am. Chem. Soc., 125, 3243 (2003).
- 5) M. Fujita, S.-Y. Yu, T. Kusukawa, H. Fumaki, K. Ogura, K. Yamaguchi: Angew. Chem., Int. Ed. Engl., 37, 2082 (1998).
- 6) S.-Y. Yu, T. Kusukawa, K. Biradha, M. Fujita: J. Am. Chem. Soc., 122, 2665 (2000).

d. Chiral Zr Catalyst Combined with Molecular Sieves

Air-stable optically active zirconium catalyst that can be stored for an extended period of time

Most catalysts easily degrade in the presence of oxygen and humidity, so must be used in inert gas systems and under strictly anhydrous conditions. These conditions have been a major obstacle to their industrialization. To solve these problems, asymmetrical Lewis acid catalysts with zirconium, as the central metal and molecular sieves (MS), a type of zeolite, were combined to develop high performance catalysts with significantly improved stability. These composite catalysts, which are air-stable and storable, are optically active Lewis acid catalysts that can be stored for an extended period of time and are easy to handle. Thus Industrialization of asymmetric catalytic reactions is highly anticipated.

[Reactions]

Asymmetric Mukaiyama Aldol Reaction

Asymmetric Hetero-Diels-Alder Reaction*

*The obtained product was treated with TFA.

	Description	Grade	Wako Cat No.	Package Size
Zr	(R)-3-I-ZrMS	for organic synthesis	097-05271	1g
	(n)-5-1-21W3	ioi organic synthesis	093-05273	5 g

[References]

- 1) S. Kobayashi, S. Saito, M. Ueno, Y. Yamashita: Chem. Commun., 2016 (2003).
- 2) S. Kobayashi, M. Ueno, S. Saito, Y. Mizuki, H. Ishitani, Y. Yamashita: Proc. Natl. Acad. Sci. USA, 101, 5476 (2004).

e. Precious Metal Powder Catalysts

The demand for high performance, high activity catalysts in the fields of specialty chemicals (fine chemicals, electronic materials, etc.) and life science (drugs, agrichemicals, and intermediates) is increasing along with tighter environmental standards and increasing requests for reduced manufacturing costs. The following 7 precious powder metal catalysts are widely suitable for reactions where these are frequently used in the above fields.

[Reactions]

Selectively hydrogenated reaction

Debenzylation reaction

Description	Purpose	Wako Cat No	Package Size	
5% Palladium on Activated Carbon, Degussa type E 101 O/W	Disproportionation, reductive	327-81641	5 g	-
(wetted with ca.55% water)	amination reaction, etc.	325-81642	25 g	_
5% Palladium on Activated Carbon, Degussa type E 101 NE/W	Deprotection, Debenzylation, etc.	324-81651	5 g	Ī,
(wetted with ca.55% water)	Deprotection, Depenzylation, etc.	322-81652	25 g	1
5% Palladium on Activated Carbon, Degussa type E 1533 R/W	Reduction of nitro group, nitrile group, etc.	321-81661	5 g	
(wetted with ca.55% water)		329-81662	25 g	
5% Palladium on Activated Carbon, Degussa type E 1002 NN/W	Hydrogenation of C=C bonds,	328-81671	5 g	
(wetted with ca.55% water)	dehalogenation reaction, etc.	326-81672	25 g	
3% Platinum on Activated Carbon, Degussa type F 1002 RCA/W	Reduction of imine, hydrogenation	325-81681	5 g	
(wetted with ca.55% water)	of heterocycles, etc.	323-81682	25 g	
1% Platinum on Activated Carbon, Degussa type F 1002 RAU/W	Selective hydrogenation (inhibits	322-81691	5 g	
(wetted with ca.55% water)	halogen desorption)	320-81692	25 g	
5% Rhodium on Activated Carbon, Degussa type G 106 NB/W	Hydrogenation of aromatic rings,	325-81701	5 g	
(wetted with ca.55% water)	reduction of nitrile, etc.	323-81702	25 g	_ (








Rh

f. Metal Scavengers (Organofunctional polysiloxanes)


Deloxan®

Precious metals dissolved in reaction liquids used in the medical industry cause serious problems such as the loss of valuable precious metals and product contamination. Deloxan® effectively solves such problems by quantitatively and efficiently adsorbing precious metals (Rh, Pd, Pt, Ir, Ru, Ag, Au) and non-precious metals (e.g. Ni, Cu, Bi, Ga, Fe) (1 ppm or less of the residual metals in the solution) as well as by allowing for the recovery of precious metals.

Deloxan® THP II Functional group: thiourea



Deloxan® MP Functional group: mercapto



Description	Purpose	Wako Cat No	Package Size
Motel Canyongov Dogusa tuno Dolovono MD	Removal of precious metals	325-83521	5 g
Metal Scavenger, Degussa type Deloxan® MP	Removal of precious metals	323-83522	25 g
Metal Scavenger, Degussa type Deloxan® THP II	Removal of precious metals	328-83511	5 g
	Removal of precious metals	326-83512	25 g

Crown ether has the property to catalyze organic reactions by forming a stable complex with metallic cation and generating a chemically stable anion. In addition, the ionic liquid can be recovered and reused because it is an organic liquid that is composed only of ions that are nonvolatile, highly conductive, and catalytically active, and thus does not blend with organic solvents. The 18-crown-6-supported ionic liquid shows the catalytic activity of the crown ether and the immiscibility with the organic solvents of ionic liquids. The use of this catalyst allows for a high yield of the objective substance in nucleophilic substitution of potassium acetate in ionic liquid solvents by benzyl bromide. To reuse the ionic liquid, the reaction liquid is extracted using liquids such as ether, and the residual ionic liquid is dried.

1-(2-Methyl-18-crown-6)-3methylimidazolium Hexafluorophosphate

[Reactions]

Acylation reaction

Fluorination reaction

$$\begin{array}{c} \text{CH}_2\text{Br} \\ \hline \\ \text{Catalyst (0.2eq)} \\ \hline \\ \text{[bmim]PF}_6 \\ \text{KF (1.5eq), 80°C / 10h} \\ \text{[or CsF (1.5eq.), rt / 12h]} \\ \end{array} \begin{array}{c} \text{CH}_2\text{F} \\ \\ \text{Yield (\%, cycle) :} \\ \text{KF : 91 (1st), 91 (2nd), 90(3rd), 91(4th), 90(5th)} \\ \text{CsF : 97 (1st), 96 (2nd), 96(3rd)} \\ \end{array}$$

<Reaction condition>

1-(2-methyl-18-crown-6)-3-methylimidazolium Hexafluorophosphate (50 mg, 0.1 mmol, 0.2 eq) and 1-butyl-3-methylimidazolium Hexafluorophosphate (6 mL) (dried by vacuum pump for 5 hours at 80 °C) are combined, to which a reaction substrate (0.5 mmol) and a reaction reagent (0.75 mmol) are added. The mixture is stirred at each temperature and time. After the reaction, the liquid is extracted using ether, then condensed to obtain the objective substance.

<Recovery and Reuse of the Catalyst>

After the reaction, dichloromethane (10 mL) and water (1 mL) are poured in the mixture (about 5 mL) of the residual 1-(2-methyl-18-crown-6)-3-methylimidazolium Hexafluorophosphate and 1-butyl-3-methylimidazolium Hexafluorophosphate obtained by ether extraction, and stirred for 4 hours at room temperature. The organic layer is left to dry for 5 hours at 80°C under reduced pressure, and is reused during for a subsequent reaction.

Description	Purpose	Wako Cat No	Package Size
1-(2-Methyl-18-crown-6)-3-methylimidazolium	for organic synthesis	134-14981	100 mg
Hexafluorophosphate		130-14983	500 mg

(R)-(-)-3,3'-Bis(4-nitrophenyl)-1,1'-binaphthyl-2,2'-diyl Hydrogen Phosphate (First paragraph is okay)

Brønsted acid-assisted chiral Lewis acid (BLA) catalyst

Enantioselective asymmetric synthesis reactions using chiral catalysts are an excellent coupling method of optically-active compounds and are thus frequently used in research. They are normally achieved using metal Lewis acid catalysts. However, some of the downfalls involve the stability of metallic complexes and the existence of residual transition metals in the compounds.

Meanwhile, asymmetric synthesis using chiral organic catalysts constitutes a research area that is gaining increased attention. In particular, chiral Bronsted acid is stable, allowing for reactions under mild conditions, as well as easy recovery of catalysts.

Asymmetric catalysts have been the exclusive organic metal catalysts until now. Chiral catalysts have rapidly become an issue of growing interest, and thus further development is anticipated in these areas.

[Reactions]

Mannich-type reaction¹⁾

Hydrophosphonylation of aldimine²³

$$Ar$$

$$O O$$

$$O O$$

$$O O$$

$$O O$$

$$Ar = 3.5 - (CF_3)_2 C_6 H_3$$

$$Catalyst (1)$$

Description	Purpose	Wako Cat No	Package Size
(R)-NP-NAP	for organic synthesis	145-07861	200 mg

[References]

- 1) T. Akiyama, J. Itoh, K. Yokota, K. Fuchibe: Angew. Chem. Int. Ed., 43, 1566-1568 (2004).
- 2) T. Akiyama, H. Morita, J. Itoh, K. Fuchibe: Org. Lett., 7, 2583-2585 (2005).
- 3) T. Akiyama, Y. Saitoh, H. Morita, K. Fuchibe: Adv. Synth. Catal., 347, in press (2005).

4. Rare Earth Triflate

Green chemistry aims to eliminate the use of auxiliary substances during organic reactions, or if unavoidable to choose substances that have little impact on the environment.

The auxiliary substances most commonly used in organic reactions are solvents.

The use of water for reactions conducted in organic solvents significantly reduces the burden on the environment. Rare earth triflate is used in various reactions as a Lewis acid catalyst that can be used in water solution.

Description	Molecular formula/CAS No.	Wako Cat No	Package Size	
Lanthanum(III) Trifluoromethanesulfonate [References] 1) S. Kobayashi: Synlett, 689 (1994). 2) J. Inanaga, Y. Yokoyama, T. Hanamoto: Tetrahedron Lett., 34. 2791 (1993). 3) J. H. Forsberg, V. T. Spaziano, T. M. Balasubramanian, G. K. Liu, S. A. Kinsley, C. Duckworth, J. J. Poteruca, P. S. Brown, J. L. Miller: J. Org. Chem., 52. 1017 (1987).		126-05231 124-05232	5 g 25 g	La
Neodymium Trifluoromethanesulfonate	(CF ₃ SO ₃) ₃ Nd=591.45 [34622-08-7]	147-08301 145-08302	5 g 25 g	Nd
Scandium(III) Trifluoromethanesulfonate [References] 1) S. Kobayashi, I. Hachiya, M. Araki, and H. Ishitani: <i>Tetrahedron Lett.</i> , 34 . 3755 (1993). 2) S. Kobayashi: <i>Synlett</i> , 689 (1994).	(CF ₃ SO ₃) ₃ Sc=492.16 [144026-79-9]	195-11391 191-11393	1 g 5 g	Sc
Thulium Trifluoromethanesulfonate	(CF ₃ SO ₃) ₃ Tm=616.14 [141478-68-4]	205-16421 203-16422	5 g 25 g	Tm
Ytterbium(III) Trifluoromethanesulfonate n-Hydrate [References] 1) S. Kobayashi, I. Hachiya, T. Takahori : Synthesis, 371 (1993). 2) S. Kobayashi, I. Hachiya : J. Org. Chem., 59. 3590 (1994). 3) S. Kobayashi : Synlett, 689 (1994).	(CF ₃ SO ₃) ₃ Yb·xH ₂ O=620.25(As an anhydride) [54761-04-5]	254-00521 252-00522	5 g 25 g	Yb

Wako Product Update No. 14

1 2

5. Maruoka Catalysts

Chiral Phase-Transfer Catalysts (PTC)

Chiral phase-transfer catalysts, which are optically active spiro ammonium salts with 2 binaphthyl rings that allow for easy molecule designing, were invented by Professor Maruoka of the Kyoto University.

(R, R)-3,4,5-Trifluorophenyl-NAS Bromide shows high catalytic activity and high enantioselectivity in asymmetric alkylation of a-amino acid derivatives¹. By using (R, R)-3,5-Bistrifluoromethylpehnyl-NAS Bromide, β -hydroxy α -amino acid derivatives, which are important chiral units of biologically active proteins, are obtained in good yield by the aldol reaction of glycine derivatives with aldehydes. It was found that erythro isomers, which are major products, were obtained with high enantioselectivity².


[Structural formula]

(R,R)-3,4,5-Trifluorophenyl-NAS Bromide Maruoka catalyst RR-Trifluorophenyl Br Form

F F Br

(S,S)-3,4,5-Trifluorophenyl-NAS Bromide

Maruoka catalyst SS -Trifluorophenyl Br Form



(R,R)-3,5-Bistrifluoromethylphenyl-NAS Bromide Maruoka catalyst RR-Bistrifluoromethylphenyl Br Form

[Reactions]

Asymmetric alkylation of α -amino acid derivative

$$\begin{array}{c} Ph \\ Ph \\ Ph \\ \end{array} \\ OBu^t + \\ Br \\ OBn \\ \end{array} \\ \begin{array}{c} R1 \\ \\ OBn \\ \end{array} \\ \begin{array}{c} (R,R)\text{-3,4,5-Trifluorophenyl-} \\ NAS \ Bromide \\ (1 \ mol\%) \\ \end{array} \\ \begin{array}{c} 1M \ citric \ acid \\ THF \\ rt, \ 10h \\ \end{array} \\ \begin{array}{c} OBu^t \\ R^1 \\ OBr \\ \end{array} \\ \begin{array}{c} OBu^t \\ R^1 \\ OBr \\ \end{array} \\ \begin{array}{c} OBu^t \\ R^1 \\ OBr \\ \end{array} \\ \begin{array}{c} OBu^t \\ R^1 \\ OBr \\ \end{array} \\ \begin{array}{c} OBu^t \\ R^1 \\ OBr \\ \end{array}$$

R¹=OBn : Yield 81%, 98% ee R¹=H : Yield 83%, 98% ee

Asymmetric aldol reaction of β -hydroxy- α -amino acid derivative

Description	Purpose	Wako Cat No	Package Size
(R, R)-3, 4, 5-Trifluorophenyl-NAS Bromide	For organic synthesis	1201-15921	100 mg
[Maruoka catalyst RR-Trifluorophenyl Br Form]		207-15923	500 mg
(S, S)-3, 4, 5-Trifluorophenyl-NAS Bromide	For organic synthesis	201-16401	100 mg
[Maruoka catalyst SS-Trifluorophenyl Br Form]		207-16403	500 mg
(R, R)-3, 5-Bistrifluoromethylphenyl-NAS Bromide	For organic synthesis	029-14921	100 mg
[Maruoka catalyst RR-Bistrifluoromethylphenyl Br Form]		025-14923	500 mg

[References]

1) T. Ooi, M. Kameda, H. Tannai, K. Maruoka: Tetrahedron Lett., 41, 8339 (2000).

Wako Pure Chemical Industries, Ltd.

http://www.wako-chem.co.jp

1-2, Doshomachi 3-Chome

Tel: 81-6-6203-3741

Fax: 81-6-6201-5964

Online Catalog:

Chuo-Ku, Osaka 540-8605, Japan

http://search.wako-chem.com

- 2) T. Ooi, M. Taniguchi, M. Kameda, K. Maruoka: Angew. Chem. Int. Ed., 41, 4542 (2002).
 - Listed products are intended for laboratory research use only, but not to be used for drug, food or human use.
 - Please visit our online catalog to search for other products from Wako; http://search.wako-chem.com
 - This brochure may contain products that cannot be exported to your country due to regulations.
 - Bulk quote requests for some products are welcomed. Please contact us.

Wako Chemicals GmbH

http://www.wako-chemicals.de

06209IBK

European Office:

Nissanstraße 2, D-41468

Neuss, Germany Tel: 49-2131-311-0 Fax: 49-2131-311100

Wako Chemicals USA, Inc.

http://www.wakousa.com
Toll-Free (U.S. only): 1-877-714-1920
Head Office (Richmond, VA):
Tel: 1-804-714-1920/ Fax: 1-804-271-7791
Los Angeles Sales Office (Irvine, CA):
Tel: 1-949-679-1700 / Fax: 1-949-679-1701
Boston Sales Office (Cambride, MA):

Tel: 1-617-354-6772 / Fax: 1-617-354-6774