

Code No.01-0152

SPECIFICATION

Alumina, Activated for Column Chromatography

| REQUIREMENT | SPECIFICATION |
|---------------------------------|---------------|
| Appearance | White powder |
| Particle size (pass 75 μ m) | min.90% |
| Precipitation volume | 1.1~1.7mL/g |
| Hexane-soluble matter | max.0.01% |
| pH(100g/L aqueous slurry, 25°C) | 9~11 |
| Mercury (Hg) | max.0.1ppm |
| Adsorption power | to pass test |