

Code No.01-0449

## SPECIFICATION

## Antimony(III) Chloride

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay Solubility in dil. hydrochloric acid Sulfur compounds (as SO <sub>4</sub> ) Sodium (Na) Copper (Cu) Lead (Pb) Arsenic (As) Iron (Fe)	White ~ slightly yellow, crystals ~ crystalline powder min.98.0%  to pass test  max.0.01%  max.0.005%  max.5ppm  max.0.003%  max.0.002%  max.0.002%

JIS K 8400-1994(\*:Additional test performed by Wako but not required by JIS)