

Code No.01-1122

SPECIFICATION

Diammonium Adipate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 0.1mol/l solution at 25°C	6.5~7.5
Residue after ignition (as sulfate)	max.0.05%
Chloride (Cl)	max.0.001%
Nitrate	to pass test (NO ₃ :max.about 0.003%)
Sulfate (SO ₄)	max.0.001%
Iron (Fe)	max.5ppm
Assay	min.98.0%