

Code No.01-1541

SPECIFICATION

L(+)-Ascorbic Acid 2-Sulfate Disodium Salt Dihydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, crystalline powder ~ powder to pass test
Solubility in water	
Specific rotation [α]D20(c=1,H2O) (after drying)	+71~+75°
Loss on drying at 105°C	max.5.0%
pH of a 100g/l solution at 25°C	6.0~8.0
Assay (after drying)	min.95.0%