

Code No.01-1581

## SPECIFICATION

## Agar

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Loss on drying at 105°C Residue after ignition (as sulfate) Lead (Pb) Iron (Fe) Gelatin Iodine-consuming substances Jelly strength	White ~ yellowish gray, crystalline powder ~ powder 14~23%(mass/mass) max.4.0%(mass/mass) max.0.005%(mass/mass) max.0.01%(mass/mass) to pass test to pass test 400~800g/cm <sup>2</sup>

JIS K 8263:2020(\*:Additional test performed by Wako but not required by JIS)This product is powder.