

# SPECIFICATION

## Acetone

for Spectrochemical Analysis

| REQUIREMENT            | SPECIFICATION          |
|------------------------|------------------------|
| Appearance             | Colorless clear liquid |
| Density (20°C)         | 0.789~0.792g/mL        |
| Refractive index 20°C  | 1.358~1.360            |
| Absorbance (330nm)     | max.0.80               |
| Absorbance (340nm)     | max.0.07               |
| Absorbance (350~400nm) | max.0.01               |
| Water                  | max.0.2%               |
| Fluorescence test      | to pass test           |
| Assay (capillary GC)   | min.99.9%              |