

SPECIFICATION

1 mol/L Ammonium Formate Solution

for HPLC

| REQUIREMENT | SPECIFICATION |
|------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Absorbance (240nm) | max.0.40 |
| Absorbance (254nm) | max.0.03 |
| Absorbance (270~400nm) | max.0.01 |
| pH at 25°C | 6.3~6.7 |
| Fluorescence test | to pass test |
| Concentration (20°C) | 0.95~1.05mol/L |