

Code No.01-2154

SPECIFICATION

Albumin, Human, recombinant expressed in plants
for Cell Culture

REQUIREMENT	SPECIFICATION
Appearance	White ~ brown, crystals ~ powder
Solubility in water	to pass test
Loss on drying over silica gel in vacuum	under investigation
pH of a 100g/L solution at 25°C	6.5~8.5
Bacterial endotoxins	under investigation
Total nitrogen (N)	under investigation
Assay (SDS-PAGE)	under investigation

Source:Plant expressed human serum albumin