

Code No.01-2257

## SPECIFICATION

## (E)-Asarone

for the Japanese Pharmacopoeia General Tests

REQUIREMENT	SPECIFICATION
Description (JP method)	White powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
* Loss on drying over silica gel in vacuum	under investigation
* Assay (HPLC)	min.98.0%

This product is suited for (E)-Asarone of general tests on The Japanese Pharmacopoeia.(\*:Additional test performed by Wako)