

## SPECIFICATION

## 5 Anion Mixture Standard Solution

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Fluoride ion (F-) (20°C)	95.0~105.0mg/L
Chloride ion (Cl-) (20°C)	190.0~210.0mg/L
Nitrite nitrogen (as N) (20°C)	9.5~10.5mg/L
Chlorate ion (ClO <sub>3</sub> <sup>-</sup> ) (20°C)	95.0~105.0mg/L
Nitrate nitrogen (as N) (20°C)	95.0~105.0mg/L