

SPECIFICATION

Amino Acids Mixture Standard Solution, Type H (High Range)
for Amino Acid Automated Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
L-Aspartic Acid	2.38~2.62 μ mol/mL
L-Threonine	2.38~2.62 μ mol/mL
L-Serine	2.38~2.62 μ mol/mL
L-Glutamic Acid	2.38~2.62 μ mol/mL
Glycine	2.38~2.62 μ mol/mL
L-Alanine	2.38~2.62 μ mol/mL
L-Valine	2.38~2.62 μ mol/mL
L-Cystine	2.38~2.62 μ mol/mL
L-Methionine	2.38~2.62 μ mol/mL
L-Isoleucine	2.38~2.62 μ mol/mL
L-Leucine	2.38~2.62 μ mol/mL
L-Tyrosine	2.38~2.62 μ mol/mL
L-Phenylalanine	2.38~2.62 μ mol/mL
Ammonium Chloride	2.38~2.62 μ mol/mL
L-Lysine	2.38~2.62 μ mol/mL
L-Histidine	2.38~2.62 μ mol/mL
L-Arginine	2.38~2.62 μ mol/mL
L-Proline	2.38~2.62 μ mol/mL