

Code No.01-2838

SPECIFICATION

Albumin, from Bovine Serum, Fatty Acid/Protease Free, Low IgG for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	Slightly greenish yellow ~ brown, crystals ~ crystalline powder to pass test
Solubility in water	
Loss on drying at 135°C	max.5.0%
pH of a 100g/L solution at 25°C	6.5~7.5
IgG	under investigation
Fatty acid	to pass test
Protease	to pass test
Assay (electrophoresis)	min.98%