

SPECIFICATION

L-Ascorbic Acid Phosphate Magnesium Salt n-Hydrate

| REQUIREMENT | SPECIFICATION |
|--|---|
| Appearance | White ~ slightly yellow, powder or mass |
| Identification | to pass test |
| Solubility in water | to pass test |
| Loss on drying | max.29% |
| Free phosphoric acid | max.1% |
| Heavy metals | max.20ppm |
| Arsenic (as As ₂ O ₃) | max.4ppm |
| Assay (HPLC) | min.90.0% |
| Bacterial endotoxins | less than 10EU/g |