

Code No.02-0219

SPECIFICATION

Boric Acid

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (H ₃ BO ₃)	min.99.5%(mass/mass)
Solubility in water	to pass test
Solubility in ethanol	to pass test
Residue after ignition with methanol (as sulfate)	max.0.02%(mass/mass)
Chloride (Cl)	max.2ppm(mass/mass)
Phosphate (PO ₄)	max.5ppm(mass/mass)
Sulfate (SO ₄)	max.0.001%(mass/mass)
Calcium (Ca)	max.0.001%(mass/mass)
Lead (Pb)	max.5ppm(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.2ppm(mass/mass)

JIS K 8863:2022(*:Additional test performed by Wako but not required by JIS)