

SPECIFICATION

Butyl Acetate
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (CH ₃ COOCH ₂ CH ₂ CH ₂ CH ₃) (GC) Density (20°C) Refractive index 20°C Water Residue after evaporation Acidity (as CH ₃ COOH) 1-Butanol [CH ₃ (CH ₂) ₃ OH] (GC) Butyl formate [HCOO(CH ₂) ₃ CH ₃] (GC) 2-Methylpropyl acetate [CH ₃ COOCH ₂ (CH ₃) ₂] (GC) Substances darkened by sulfuric acid	Colorless clear liquid min.99.0%(mass/mass) 0.878~0.883g/ml 1.392~1.397 max.0.1%(mass/mass) max.0.002%(mass/mass) max.0.003%(mass/mass) max.0.3%(mass/mass) max.0.1%(mass/mass) max.0.5%(mass/mass) to pass test

JIS K 8377:2014

(*:Additional test performed by Wako but not required by JIS)