

SPECIFICATION

Barium Perchlorate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [Ba(ClO ₄) ₂] Solubility in water pH of a 50g/L solution at 25°C Chloride (Cl) Nitrate (NO ₃) Calcium (Ca) Strontium (Sr) Lead (Pb) Iron (Fe)	White, crystals ~ crystalline powder min.95.0%(mass/mass) to pass test 5.0~8.0 max.0.005%(mass/mass) max.0.003%(mass/mass) max.0.01%(mass/mass) max.0.1%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass)

JIS K 9551 : 2020(*:Additional test performed by Wako but not required by JIS)