

SPECIFICATION

Bis[(+)-tartrato]diantimonate(III) Dipotassium Trihydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (C ₈ H ₄ K ₂ O ₁₂ Sb ₂ ·3H ₂ O) Solubility in water pH of a 50g/l solution at 25°C Chloride (Cl) Sulfate (SO ₄) Copper (Cu) Calcium (Ca) Zinc (Zn) Lead (Pb) Arsenic (As) Iron (Fe)	White, crystals ~ crystalline powder min.99.5%(mass/mass) to pass test 3.8~4.5 max.0.001%(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass) max.0.005%(mass/mass) max.0.001%(mass/mass) max.0.005%(mass/mass) max.5ppm(mass/mass) max.5ppm(mass/mass)

JIS K 8533:2012

(*:Additional test performed by Wako but not required by JIS)