

SPECIFICATION

Butyl Acetate, Dehydrated
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.878~0.883g/mL
Refractive index 20°C	1.392~1.397
Water	max.0.005%
Acidity (as CH ₃ COOH)	max.0.003%
1-Butanol	max.0.3%
Butyl formate	max.0.1%
Isobutyl acetate	max.0.5%
Assay (capillary GC)	min.99.0%