

SPECIFICATION

Baicalin

for the Japanese Pharmacopoeia Crude Drugs Test (for TLC)

REQUIREMENT	SPECIFICATION
Description (JP method)	Light yellow, crystals or powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
* Solubility in methanol	to pass test
* Specific absorbance E1%1cm (277nm) (after drying)	589~639
* Water	measured value
* Assay (HPLC)	min.98.0%

*:Additional test performed by Wako