

Code No.03-0038

## SPECIFICATION

## Calcium Carbonate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White powder
Assay (CaCO <sub>3</sub> ) (after drying)	min.99.5%(mass/mass)
Solubility in hydrochloric acid	to pass test
Loss on drying at 105°C	max.0.3%(mass/mass)
Chloride (Cl)	max.0.005%(mass/mass)
Nitrate (NO <sub>3</sub> )	to pass test
Phosphate (PO <sub>4</sub> )	max.0.002%(mass/mass)
Sulfate (SO <sub>4</sub> )	max.0.005%(mass/mass)
Sodium (Na)	max.0.01%(mass/mass)
Potassium (K)	max.0.005%(mass/mass)
Magnesium (Mg)	max.0.01%(mass/mass)
Strontium (Sr)	max.0.05%(mass/mass)
Barium (Ba)	max.0.005%(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Ammonium (NH <sub>4</sub> )	max.0.05%(mass/mass)

JIS K 8617:2007(\*:Additional test performed by Wako but not required by JIS)