

Code No.03-0174

## SPECIFICATION

## Diammonium Cerium(IV) Nitrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Yellowish red, crystals ~ crystalline powder
Assay	min.95.0%
Solubility in water	to pass test
Chloride (Cl)	max.0.001%
Sulfate (SO <sub>4</sub> )	max.0.005%
Sodium (Na)	max.0.05%
Potassium (K)	max.0.05%
Copper (Cu)	max.0.002%
Lead (Pb)	max.0.002%
Iron (Fe)	max.0.002%

JIS K 8556-1995

(\*:Additional test performed by Wako but not required by JIS)