

Code No.03-0218

## SPECIFICATION

## Sodium p-Toluenesulfonchloramide Trihydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly yellow, crystals ~ crystalline powder
Assay (C <sub>7</sub> H <sub>7</sub> ClNaO <sub>2</sub> S · 3H <sub>2</sub> O)	98.0~104.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	8.0~11.0
Substances darkened by sulfuric acid	to pass test

JIS K 8318:2023(\*:Additional test performed by Wako but not required by JIS)