

SPECIFICATION

Chloroform

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (CHCl ₃) (GC) Density (20°C) Water Residue after evaporation Acidity (as HCl) Chloride (Cl) Free chlorine (as Cl) Carbonyl compounds (as CO) Substances darkened by sulfuric acid Suitability for use with dithizone Ethanol (C ₂ H ₅ OH) (GC)	Colorless clear liquid min.99.0%(area/area) 1.471~1.484g/mL max.0.05%(mass/mass) max.0.001%(mass/mass) max.0.001%(mass/mass) max.1ppm(mass/mass) max.1ppm(mass/mass) max.0.005%(mass/mass) to pass test to pass test 0.3~1.0%(mass/mass)

JIS K 8322:2020(*:Additional test performed by Wako but not required by JIS)