

SPECIFICATION

Copper(I) Chloride

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (CuCl) Solubility in acid Sulfate (SO ₄) Nitrogen compounds (as N) Sodium (Na) Potassium (K) Copper(II) (Cu ²⁺) Calcium (Ca) Lead (Pb) Arsenic (As) Iron (Fe)	White ~ slightly green, crystalline powder min.95.0%(mass/mass) to pass test max.0.2%(mass/mass) max.0.005%(mass/mass) max.0.05%(mass/mass) max.0.01%(mass/mass) max.0.5%(mass/mass) max.0.01%(mass/mass) max.0.02%(mass/mass) max.1ppm(mass/mass) max.0.005%(mass/mass)

JIS K 8138:2013(Amendment 1)

(*:Additional test performed by Wako but not required by JIS)