

Code No.03-0490

SPECIFICATION

Cupferron

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|-------------------------------------|--|
| * Appearance | White ~ pale yellowish brown, crystals ~ |
| Solubility in water | crystalline powder |
| Specific absorbance[E1%1cm(281nm)] | to pass test |
| Loss on drying | min.515 |
| Residue after ignition (as sulfate) | max.1.0%(mass/mass) |
| Sensitivity | max.0.05%(mass/mass) |
| | to pass test |

JIS K 8289:2012(*:Additional test performed by Wako but not required by JIS)