

Code No.03-0490

## SPECIFICATION

## Cupferron

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Solubility in water Specific absorbance[E1%1cm(281nm)] Loss on drying Residue after ignition (as sulfate) Sensitivity	White ~ pale yellowish brown, crystals ~ crystalline powder to pass test  min.515  max.1.0%(mass/mass)  max.0.05%(mass/mass)  to pass test

JIS K 8289:2012(\*:Additional test performed by Wako but not required by JIS)