

Code No.03-0490

SPECIFICATION

Cupferron

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|---|--|
| * Appearance Solubility in water Specific absorbance[E1%1cm(281nm)] Loss on drying Residue after ignition (as sulfate) Sensitivity | White ~ pale yellowish brown, crystals ~ crystalline powder to pass test min.515 max.1.0%(mass/mass) max.0.05%(mass/mass) to pass test |

JIS K 8289:2012(*:Additional test performed by Wako but not required by JIS)