

Code No.03-0506

## SPECIFICATION

## Cyclohexanol

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (C <sub>6</sub> H <sub>12</sub> O) (GC)	min.98.0%(mass/mass)
Freezing point	min.20°C
Density (25°C)	0.943~0.953g/ml
Water	max.0.1%(mass/mass)
Acidity (as C <sub>6</sub> H <sub>5</sub> COOH)	max.0.03%(mass/mass)
Cyclohexanone (C <sub>6</sub> H <sub>10</sub> O) (GC)	max.0.5%(mass/mass)

JIS K 8462:2011

(\*:Additional test performed by Wako but not required by JIS)