

Code No.03-0598

SPECIFICATION

Citral (mixture of cis- and trans-)

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellow, clear liquid
Solubility in ethanol	to pass test
Density (20°C)	0.884~0.892g/ml
Refractive index 20°C	1.486~1.490
Acidity (as C ₁₀ H ₁₆ O ₂)	max.0.5%
Assay (capillary GC) (cis- + trans-)	min.95.0%