

Code No.03-1350

SPECIFICATION

Calcium Propionate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, granules ~ powder
Solubility in dil.acetic acid	to pass test
Loss on drying at 120°C	max.3.0%
Acidity (as C ₂ H ₅ COOH)	max.0.1%
Alkalinity [as Ca(OH) ₂]	max.0.2%
Heavy metals (as Pb)	max.0.001%
Arsenic (As)	max.2ppm
Assay (after drying at 120°C)	min.98.0%