

Code No.03-1416

SPECIFICATION

Calcium (+)-Pantothenate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White powder
Solubility in water	to pass test
Specific rotation [α]D20(c=5,H2O)	+26.0~+28.0°
Loss on drying at 105°C	max.5.0%
Chloride (Cl)	max.0.02%
Sulfate (SO4)	max.0.03%
Heavy metals (as Pb)	max.0.002%
Assay (after drying)	min.98.0%