

Code No.03-2033

SPECIFICATION

(E)-Chlorogenic Acid

for the Japanese Pharmacopoeia Crude Drugs Test (for Thin-Layer Chromatography)

| REQUIREMENT | SPECIFICATION |
|--|---------------------|
| Description (JP method) | White powder |
| Purity (JP method) | - |
| Related substances (JP method) | to pass test |
| * Solubility in methanol | to pass test |
| * Specific rotation [α]D20(c=0.4,H2O) | under investigation |
| * Water | under investigation |
| * Assay (HPLC) | min.98.0% |

*:Additional test performed by Wako