

## SPECIFICATION

Chloroform, Super Dehydrated, Amylene added  
for Organic Synthesis

| REQUIREMENT                          | SPECIFICATION                    |
|--------------------------------------|----------------------------------|
| Appearance                           | Colorless clear liquid           |
| Water                                | max.0.001%                       |
| Residue after evaporation            | max.0.001%                       |
| Acidity (as HCl)                     | max.0.001%                       |
| Chloride (Cl)                        | max.1ppm                         |
| Free chlorine                        | to pass test (Cl:max.about 1ppm) |
| Substances darkened by sulfuric acid | to pass test                     |
| Assay (capillary GC)                 | min.99.5%                        |
| Amylene (stabilizer)                 | 100~250ppm                       |