

Code No.03-2540

## SPECIFICATION

### Cap B Solution [Acetic Anhydride-2,6-Lutidine-Acetonitrile (2:3:5)] for Nucleic Acid Synthesis

| REQUIREMENT           | SPECIFICATION   |
|-----------------------|-----------------|
| Appearance            | to pass test    |
| Density (20°C)        | 0.883~0.889g/mL |
| Refractive index 20°C | 1.398~1.402     |
| Impurities            | to pass test    |
| Acetonitrile          | 49.0~51.0vol%   |
| Acetic anhydride      | 19.0~21.0vol%   |
| 2,6-Lutidine          | 29.0~31.0vol%   |