

Code No.04-0276

## SPECIFICATION

## N,N-Dimethylaniline

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellow, clear liquid
Density (20°C)	0.954~0.958g/mL
Refractive index 20°C	1.556~1.560
Water	max.0.2%
Aniline (C <sub>6</sub> H <sub>5</sub> NH <sub>2</sub> )	max.0.03%
N-Methylaniline (C <sub>6</sub> H <sub>5</sub> NHCH <sub>3</sub> )	max.0.3%
Assay (capillary GC)	min.99.0%