

Code No.04-0392

## SPECIFICATION

## 1,5-Diphenylcarbonohydrazide

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly red, crystalline powder ~ powder
Solubility in acetone	to pass test
Melting point	171~175°C
Loss on drying over silica gel in vacuum	max.1.0%(mass/mass)
Residue after ignition (as sulfate)	max.0.05%(mass/mass)
Suitability for chromium determination	-
Sample test	min.0.42
Blank test	max.0.01

JIS K 8488:2011

(\*:Additional test performed by Wako but not required by JIS)